

■ 显微外科器械

Microsurgery Instruments

■ 基础外科器械

General Instruments



1

Tel: +86 517 3160500

Fax: +86 517 3160155

E-mail: info@tisurg.com

Table of Contents

目录

■ Microsurgery /General Instruments 显微外科/基础外科器械

■ Forceps 镊子	
Micro Forceps(Tying、Toothed、Ring) 显微镊	2-3
Adson Tissue Forceps Adson镊	3-4
Tissue Forceps(Bishop-Harmon、Tying、Graefe、Thumb) 组织镊	4
Jewelers Forceps 精细镊	5
■ Needle Holders 持针钳	
Micro Needle Holder 显微持针钳	6-8
Finger Ring Style Needle Holder 环形指圈持针钳	9
Olsen-Hegar Needle Holder 带剪持针钳	9
■ Scissors 剪刀	
Micro Scissors(Micro、Vannas、Noyes、Tenotomy、Jacobson) 显微剪	10
Super Cut Scissors 快口剪刀	11
Finger Ring Style Scissors(Metzenbaum、Mayo) 环形指圈剪刀	11
■ Haemostatic Forceps 止血钳	
Mosquito Forceps 蚊式钳	12
Haemostatic Forceps(Crile、Leriche、Bengolea、Kelly、O'shaughnessy) 止血钳	13-14

- Microsurgery Instruments
- General Instruments

显微外科器械
基础外科器械

Micro Forceps 显微镊

- T1-1101 Straight, 0.1mm tips
- T1-1102 Straight, 0.2mm tips
- T1-1103 Straight, 0.4mm tips



Tungsten carbide coated platform, round handle, 15cm

Micro Forceps 显微镊

- T1-1111 Straight, 0.2mm tips
- T1-1112 Straight, 0.3mm tips
- T1-1113 Straight, 0.6mm tips

- T1-1121 Curved, 0.2mm tips
- T1-1122 Curved, 0.3mm tips
- T1-1123 Curved, 0.6mm tips



Tungsten carbide coated platform, round handle, 15cm

Micro Forceps 显微镊

- | | |
|------------------------------------|--------------------------------------|
| T1-1131 Straight, 0.2mm tips, 18cm | T1-1141 Curved, 0.2mm tips, 18cm |
| T1-1132 Straight, 0.3mm tips, 18cm | T1-1142 Curved, 0.3mm tips, 18cm |
| T1-1133 Straight, 0.7mm tips, 18cm | T1-1143 Curved, 0.7mm tips, 18cm |
| T1-1151 Straight, 0.3mm tips, 20cm | T1-1161 Straight, 0.3mm tips, 23cm |
| T1-1152 Curved, 0.3mm tips, 20cm | T1-1162 Curved, 0.3mm tips, 23cm |
| T1-1153 Straight, 0.7mm tips, 20cm | T1-1163 Straight, 0.7mm tips, 23cm |
| T1-1154 Curved, 0.7mm tips, 20cm | T1-1164 Curved, 0.7mm tips, 23cm |
| T1-1171 Straight, 0.5mm tips, 25cm | T1-1181 Straight, 0.7mm tips, 27.5cm |
| T1-1172 Curved, 0.5mm tips, 25cm | T1-1182 Curved, 0.7mm tips, 27.5cm |
| T1-1173 Straight, 1.0mm tips, 25cm | T1-1191 Straight, 1.0mm tips, 30cm |
| T1-1174 Curved, 1.0mm tips, 25cm | T1-1192 Curved, 1.0mm tips, 30cm |

T1-1132



T1-1142



Tungsten carbide coated platform, round handle

Micro Forceps 显微镊

- T1-1155 Straight, 20cm
- T1-1165 Straight, 23cm
- T1-1175 Straight, 25cm
- T1-1185 Straight, 27.5cm

T1-1185




10mm tungsten carbide coated platform, round handle

Micro Toothed Forceps 带齿显微镊

- T1-1201 With platform, straight, 0.3mm teeth, 15cm
- T1-1202 With platform, curved, 0.3mm teeth, 15cm
- T1-1203 With platform, straight, 0.5mm teeth, 15cm
- T1-1204 With platform, curved, 0.5mm teeth, 15cm
- T1-1205 Without platform, straight, 0.5mm teeth, 15cm
- T1-1206 Without platform, curved, 0.5mm teeth, 15cm

- T1-1231 With platform, straight, 0.3mm teeth, 18cm
- T1-1232 With platform, curved, 0.3mm teeth, 18cm
- T1-1233 With platform, straight, 0.5mm teeth, 18cm
- T1-1234 With platform, curved, 0.5mm teeth, 18cm
- T1-1235 Without platform, straight, 0.5mm teeth, 18cm
- T1-1236 Without platform, curved, 0.5mm teeth, 18cm



T1-1233



Tungsten carbide coated platform, 1X2teeth, round handle

Micro Ring Tip Forceps 环形头显微镊

- | | |
|----------------------------------|------------------------------------|
| T1-1301 1.0mm ring tips, 15cm | T1-1331 1.0mm ring tips, 18cm |
| T1-1302 1.3mm ring tips, 15cm | T1-1332 1.3mm ring tips, 18cm |
| T1-1351 1.0mm ring tips, 20cm | T1-1361 1.0mm ring tips, 23cm |
| T1-1352 1.3mm ring tips, 20cm | T1-1362 1.3mm ring tips, 23cm |
| T1-1353 1.5mm ring tips, 20cm | T1-1363 1.5mm ring tips, 23cm |
| T1-1371 1.0mm ring tips, 25cm | T1-1381 1.0mm ring tips, 27.5cm |
| T1-1372 1.3mm ring tips, 25cm | T1-1382 1.3mm ring tips, 27.5cm |
| T1-1373 1.5mm ring tips, 25cm | T1-1383 1.5mm ring tips, 27.5cm |

T1-1353



Tungsten carbide coated tips, round handle

Adson Tissue Forceps Adson 组织镊

- | | |
|---|---|
| T1-1401 Without platform, straight, 12cm | T1-1406 Without platform, straight, 15cm |
| T1-1402 With platform, straight, 12cm | T1-1407 With platform, straight, 15cm |
| T1-1403 With platform, curved, 12cm | T1-1408 With platform, curved, 15cm |



T1-1401




Tungsten carbide coated tips, flat handle

Adson Tissue Forceps Adson 组织镊

- T1-1411** Without platform, straight, 1×2teeth, 12cm
T1-1412 With platform, straight, 1×2teeth, 12cm
T1-1413 With platform, curved, 1×2teeth, 12cm
T1-1414 Without platform, straight, 2×3teeth, 12cm

T1-1416 Without platform, straight, 1×2teeth, 15cm
T1-1417 With platform, straight, 1×2teeth, 15cm
T1-1418 With platform, curved, 1×2teeth, 15cm
T1-1419 Without platform, straight, 2×3teeth, 15cm



T1-1411



Tungsten carbide coated tips, flat handle

Adson Tissue Forceps Adson 组织镊

- T1-1421** Straight, 7×7teeth, 12cm
T1-1422 Curved, 7×7teeth, 12cm **T1-1421**
T1-1426 Straight, 7×7teeth, 15cm
T1-1427 Curved, 7×7teeth, 15cm



Bishop-Harmon Forceps 组织镊

- T1-1501** 0.5mm tips, 9cm
T1-1502 1×2teeth, 9cm



Tungsten carbide coated tips, flat handle

Tying Forceps 系线镊

- T1-1511** 0.4mm tips, with platform, 9cm
T1-1512 1×2teeth, with platform, 9cm



Tungsten carbide coated tips, flat handle


Graefe Forceps 组织镊

- T1-1521** 0.8mm serrated tips, straight
T1-1522 0.8mm 1×2teeth, straight
T1-1523 0.8mm serrated tips, curved
T1-1524 0.8mm 1×2teeth, curved

T1-1525 0.8mm serrated tips, strongly curved
T1-1526 0.8mm 1×2teeth, strongly curved
T1-1527 0.8mm serrated tips, angled
T1-1528 0.8mm 1×2teeth, angled



T1-1522



Overall length 11cm, flat handle

Thumb Forceps 组织镊

- T1-1601** 1mm tips, 15cm
T1-1602 2mm tips, 15cm
T1-1603 1×2teeth, 15cm

T1-1651 1mm tips, 20cm
T1-1652 2mm tips, 20cm
T1-1653 1×2teeth, 20cm

T1-1631 1mm tips, 18cm
T1-1632 2mm tips, 18cm
T1-1633 1×2teeth, 18cm

T1-1661 1mm tips, 24.2cm
T1-1662 2mm tips, 24.2cm
T1-1663 1×2teeth, 24.2cm



T1-1652



- Microsurgery Instruments 显微外科器械
- General Instruments 基础外科器械

Jewelers Forceps 精细镊

T1-1801 #1 Straight, 0.10x0.06mm tips, 12cm

T1-1802 #1 Straight, 0.20x0.12mm tips, 12cm



T1-1803 #2 Straight, 0.17x0.07mm tips, 12cm

T1-1804 #2 Straight, 0.34x0.14mm tips, 12cm



T1-1805 #3 Straight, 0.08x0.04mm tips, 12cm

T1-1806 #3 Straight, 0.17x0.10mm tips, 12cm

T1-1807 #3c Straight, 0.08x0.04mm tips, 11cm



T1-1808 #4 Straight, 0.06x0.02mm tips, 11cm

T1-1809 #4 Straight, 0.13x0.08mm tips, 11cm



T1-1810 #5 Straight, 0.05x0.01mm tips, 11cm

T1-1811 #5 Straight, 0.10x0.06mm tips, 11cm

T1-1812 #5 Straight, with platform, 11cm



With platform



T1-1813 #5 Angled 45°, 0.05x0.01mm tips, 11cm

T1-1814 #5 Angled 45°, with platform, 11cm



With platform



T1-1815 #5 Angled 90°, 0.05x0.01mm tips, 11cm

T1-1816 #5 Angled 90°, with platform. 11cm



With platform

T1-1817 #7 Curved, 0.17x0.10mm tips, 11.7cm



T1-1818 0.15x0.12mm tips, 15.5cm

T1-1819 With platform, 15.5cm



With platform

