

# 46

## Lipo Suctions Instruments

Cannulas  
Micro Fat Transfer  
Decent Stands  
Sterilization Box  
LipoSuction Gun

Edition  
**A17**



# HANDLES / CANNULA

**1****SIMPLE LUER LOCK**

*"Most Economical and commonly used cannula, Simply Stainless Steel material."*

**2****SUPER LUER LOCK**

*"Helpful for distribute fluid, specifically for the infusion of Tumescent solutions and for hydraulic dissections"*

**3****TOMMEY**

*"Helpful for distribute fluid, specifically for the infusion of Tumescent solutions and for hydraulic dissections"*

**4****HANDLE FIX**

*"Helpful for distribute fluid, specifically for the infusion of Tumescent solutions and for hydraulic dissections"*

**5****S. S Golden Handle**

*"Helpful for distribute fluid, specifically for the infusion of Tumescent solutions and for hydraulic dissections"*



### INFUSION



#### 1 FLUID INFUSION CANNULA

*"Helpful for distribute fluid, specifically for the infusion of Tumescent solutions and for hydraulic dissections"*

#### Available Sizes

Dia	Length	Holes
1.5mm		08
2.0mm		08
2.5mm		08
3.0mm		08
<b>Recommendably lengths are 12cm, 15cm, 20cm, 25cm, 30cm</b>		

### INJECTORS



#### 1 ROUND Blunt WIDE Injector

*"Specifically used for soft Tissue and can handle small volume injector"*

#### Available Sizes

Dia	Length	Holes
0.7mm		01
0.9mm		01
1.2mm		01
1.2mm		01 Curved
1.4mm		01 Curved
1.6mm		01 Curved
<b>Recommendably lengths are 4cm, 5cm, 7cm, 9cm, 12cm, 15cm, 20cm</b>		



#### 2 BLUNT INJECTOR

*"Precisely for Micro fat Injector, Useful for periorbital region and vermilion skin boarder for the lips"*

0.7mm		01
0.9mm		01
1.2mm		01
1.6mm		01
<b>Recommendably lengths are 4cm, 5cm, 7cm, 9cm</b>		



#### 3 Spoon Injectors

*"Precisely for large fat Injectors, Useful for breast and buttock workhorse"*

1.2mm		01
1.4mm		01
1.6mm		01
2.0mm		01
2.4mm		01
3.0mm		01
1.4mm		01 Curved
1.6mm		01 Curved



#### 4 SUPER Facial Injectors

*"Precisely for Intradermal injection of Micro fat."*

<b>Recommendably lengths are 7cm, 10cm, 15cm, 20cm, 25cm</b>		
0.7mm		05
09mm		05
2.4mm		05
3.0mm		05
3.7mm		05






## INJECTOR Cannulas LIPOSUCTION










	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	S.S Golden Handle
	46.105115	46.105215	46.105315	46.105415	46.105515
	46.105120	46.205120	46.305120	46.405120	46.105520
	46.105125	46.205125	46.305125	46.405125	46.105525
	45.105130	46.205130	46.305130	46.405130	46.10530

	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	S.S Golden Handle
	46.111107	46.111207	46.111307	46.111407	46.111507
	46.111109	46.111209	46.111309	46.111409	46.111509
	46.111112	46.111212	46.111312	46.111412	46.111512
	46.112112	46.112212	46.112312		
	46.112114	46.112214	46.112314		
	46.112116	46.112216	46.112316		
	46.113107	46.113207	46.113307	46.113407	46.113507
	46.113109	46.113209	46.113309	46.113409	46.113509
	46.113112	46.113213	46.113313	46.113413	46.113513
	46.113116	46.113216	46.113316	46.113416	46.113516
		46.115214			
	46.115114	46.115214			
	46.115116	46.115216			
	46.115120	46.115220			
	46.115124	46.115224	46.115324		46.115524
	46.115130	46.115230	46.115330		46.115524
	46.116114	46.116214			
	46.116116	46.116216			
	46.117107				
	46.117109				
			46.117324		46.117524
			46.117330		46.117530
			46.117337		46.117537



# HARVESTER

### Available Sizes


	Dia	Length	Holes
<b>1</b>  <b>BEVELED EDGE HARVESTER</b> <i>"Six Beveled Edges, It deliver maximum Tissue volume in minimal time"</i>	2.1mm		05
	<b>Recommendably lengths are 10cm, 20cm</b>		
<b>2</b>  <b>MICRO HARVETER</b> <i>"Precisely for Selective and safe removal, for SVF-Rich (Stromal Vascular Fraction) fat Layer in deep dermis"</i>	2.1mm		04
	<b>Recommendably lengths are 10cm, 20cm</b>		
<b>3</b>  <b>THREE PORT HARVESTER</b> <i>"Usually used for Multi plane, adipose tissue harvesting, very popular for small volume procedures"</i>	1.6mm		03
	2.0mm		03
	2.5mm		03
	3.0mm		03
	3.5mm		03
	4.0mm		03
	5.0mm		03
	<b>Recommendably lengths are 10cm, 12cm, 15cm, 20cm, 25cm, 30cm</b>		
<b>4</b>  <b>LITTLE Harvester</b> <i>"It's a less aggressive harvester, 6 port configuration"</i>	2.1mm		06
	<b>Recommendably lengths are 10cm, 12cm, 15cm</b>		
<b>5</b>  <b>QUICK Harvester</b> <i>"Six port, less aggressive, and maximum tissue volume in minimal time"</i>	2.5mm		06
	3.0mm		06
	3.7mm		06
<b>Recommendably lengths are 15cm, 20cm, 25cm, 30cm</b>			
<b>6</b>  <b>30 Holes Harvester</b> <i>Quick harvesting, for small area, and also help full for reducing blocking</i>	2.5cm		30
	3.0cm		30
<b>Recommendably lengths are 15cm, 20cm, 25cm</b>			
<b>7</b>  <b>MICROPOR SPEED Harvester</b> <i>"Ports are beveled and sharp elevated positioned symmetrically to create precise Impact"</i>	2.5mm		10
	3.0mm		10
	3.7mm		10
<b>Recommendably lengths are 15cm, 20cm, 25cm</b>			



**HARVESTER Cannula**  
LIPOSUCTION








	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	S.S Golden Handle
	46.131121	46.131221			46.131521
	46.133121	46.133221			46.133521
	46.135116	46.135216			
	46.135120	46.135220			
	46.135125	46.135225	46.135325		46.135525
	46.135130	46.135230	46.135330		46.135530
			46.115335		46.135535
			46.115340		46.135540
			46.115350		46.135550
	46.137107	46.137207			46.137507
	46.139125	46.139225	46.139325	46.139425	46.139525
	46.139130	46.139230	46.139330	46.139430	46.139530
	46.139137	46.139237	46.139337	46.139437	46.139537
	46.141125	46.141225	46.141325		
	46.141130	46.141230	46.141330	46.141430	46.141530
	46.143125	46.143225	46.143325	46.143425	
	46.143130	46.143230	46.143330	46.143430	46.143530
				46.143437	46.143537



# HARVESTER

### Available Sizes


	Dia	Length	Holes
<b>8</b>  <b>MULTIPOINT HARVESTER</b> <i>"Ports are beveled, 12 to 20 port, efficient harvester both forward and backward movement"</i>	2.4mm		12-20
	3.0mm		12-20
	3.7mm		12-20
	<b>Recommendably lengths are 10cm, 20cm</b>		
<b>9</b>  <b>MULTIPOINT AGGRESSIVE Harvester</b> <i>"Multiport Harvest, with elevated sharp micro port, effective for fibrous tissue"</i>	2.0mm		10
	2.4mm		10
	3.0mm		10
	3.7mm		10
	<b>Recommendably lengths are 15cm, 20cm, 25cm, 30cm, 35cm</b>		
<b>10</b>  <b>TRIPLE PORT Harvester</b> <i>"Three opening ports in a linear pattern, it is Best useful for adipose tissues,"</i>	2.0mm		04
	2.4mm		04
	3.0mm		04
	3.7mm		04
	4.0mm		04
	4.6mm		04
	<b>Recommendably lengths are 15cm, 20cm, 25cm, 30cm, 35cm</b>		
<b>11</b>  <b>QUAD PORT Harvester</b> <i>"Four openings in a linear pattern. Distal end of each opening is raised for etching of subcutaneous tissue"</i>	3.0mm		04
	<b>Recommendably lengths are 20cm</b>		
<b>12</b>  <b>LANTERN Shaped Cannula</b> <i>"Lantern Shaped cannula, with three prot, useful for Fibrous areas, skin irregularities, as an equalizer, injector and popular in Gynecomastia procedures"</i>	3.0mm		03
	3.7mm		03
	4.0mm		03
	5.0mm		03
	<b>Recommendably lengths are 9cm, 25cm, 30cm</b>		



# HARVESTER Cannula LIPOSUCTION








	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	Adopter Luer Lock
	46.145124	46.145224	46.145324	46.145424	46.145524
	46.145130	46.145230	46.145330	46.145430	46.145530
				46.145437	46.145537
	46.147120	46.147220	46.147320	46.147420	
	46.147124	46.147224	46.147324	46.147424	46.147524
	46.147130	46.146230	46.146330	46.147430	46.147530
				46.146437	46.147537
	46.149120	46.149220			
	46.149124	46.149224	46.149324	46.149424	
	46.149130	46.149230	46.149330	46.149430	46.149530
	46.149137	46.149237	46.149337	46.149437	46.149537
			46.149340	46.149440	46.149540
			46.149346	46.149446	46.149546
				46.151430	46.151530
		46.153230		46.153430	
				46.153437	46.153537
			46.153340	46.153440	46.153540
				46.153450	46.153550





# CANNULAS

### Available Sizes


	Dia	Length	Holes
<b>1</b>  <b>COBBRA Lipo Cannula</b> <i>"Tapered parallel two ports, working in heavy tissues- aggressive than Conventional, Pyramid and Mercedes cannulas"</i>	1.5mm		02
	2.0mm		02
	2.5mm		02
	3.0mm		02
	3.7mm		02
	4.6mm		02
	<b>Recommendably lengths are</b>		
<b>2</b>  <b>OPEN Harvester</b> <i>"Guarded fully open tip, designed specifically for tissue harvest"</i>	3.0mm		02
	<b>Recommendably lengths are</b>		
<b>3</b>  <b>MERCEDES Lip Cannula</b> <i>"Parallel three ports, for multi plane soft tissues."</i>	2.0mm		03
	2.5mm		03
	3.0mm		03
	3.7mm		03
	4.0mm		03
	4.6mm		03
	5.0mm		03
	6.0mm		03
<b>Recommendably lengths are</b>			<b>10cm, 15cm, 20cm, 25cm, 30cm, 36cm</b>
<b>4</b>  <b>ADJACENT Holes Cannulas</b> <i>"Two Parallel Hales, above a single holes, at 180 Holes, aggressive than conventional cannulas"</i>	2.0mm		03
	3.0mm		03
	4.0mm		03
	5.0mm		03
	6.0mm		03
	<b>Recommendably lengths are</b>		
<b>5</b>  <b>PYRAMID Lipo Cannula</b> <i>"Its three port cannula, at 180 degree and slightly more aggressive than conventional cannula"</i>	2.0mm		03
	2.5mm		03
	3.0mm		03
	3.7mm		03
	4.0mm		03
	4.6mm		03
	<b>Recommendably lengths are</b>		



**CANNULAS**  
LIPOSUCTION



	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	Adopter Luer Lock
		46.161215			
	46.161120	46.161220			
	46.161125	46.161225	46.161325	46.161425	46.161525
	46.161130	46.161230	46.161330	46.161430	46.161530
			46.161337	46.161437	46.161537
			46.161346		
	46.163130	46.163230	46.163330		46.163537
	46.165120	46.165220	46.165320	46.165420	
	46.165125	46.165225	46.165325	46.165425	46.165525
	46.165130	46.165230	46.165330	46.165430	46.165530
			46.165337	46.165437	46.165537
			46.165340	46.165440	46.165540
			46.165346	46.165446	46.165546
			46.165350	46.165450	46.165550
			46.165360	46.165460	46.165560
		46.167220	46.167320	46.167420	
		46.167230	46.167330	46.167430	46.167530
			46.167340	46.167440	46.167540
			46.167350	46.167450	46.167550
			46.167360	46.167460	46.167560
	46.169120	46.169220	40.169320	46.169420	
	46.169125	46.169225	46.169325	46.169425	46.169525
	46.169130	46.169230	46.169330	46.169430	46.169530
			46.169335	46.169435	46.169535
			46.169340	46.169440	46.169540
			46.169346		



## CANNULAS

6



## CONVENTIONAL Lipo Cannula

*"Deep Single tissue removal port, commonly used for tissue removal"*

7



## SPEEDO Liposculpture Cannula

*"It's an advanced Lipo-sculpture cannulas, for Quick and effective liposuctions"*

## Available Sizes


Dia	Length	Holes
0.7mm		01
0.9mm		01
1.2mm		01
1.4mm		01
1.6mm		01
2.0mm		01
2.5mm		01
3.0mm		01
3.7mm		01
4.0mm		01
4.6mm		01
<b>Recommendably lengths are 3cm, 5cm, 7cm, 9cm, 12cm, 15cm, 20cm, 25cm, 30cm, 35cm</b>		
2.5mm		
3.0mm		
3.7mm		
<b>Recommendably lengths are 15cm, 20cm, 25cm, 30cm, 35cm</b>		



**CANNULAS**  
LIPOSUCTION








	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	Adopter Luer Lock
	46.171107				
	46.171109				
	46.171112				
	46.171114				
	46.171116				
	46.171120	46.171220		46.171420	46.171520
	46.171125	46.171225	46.171325	46.171425	46.171525
	46.171130	46.171230	46.171330	46.171430	46.171530
			46.171337	46.171437	46.171537
			46.171340	46.171440	46.171540
			46.171346		
				46.173424	46.173524
	46.173130	46.173230	46.173330	46.173430	46.173530
			46.173337	46.173437	46.173537





# FILTRATOR

### Available Sizes

	Dia	Length	Holes
<p><b>1</b></p>  <p><b>Facelift Filtrators Cannula</b> <i>"It reduces bleeding and dissection time through hydraulic dissection."</i></p>	2.0mm		06
	<b>Recommendably lengths are 15cm</b>		
<p><b>2</b></p>  <p><b>FLAPPER Cannulas</b> <i>"It is used undermining dissector, followed by harvesting and/or injection"</i></p>	2.0mm		01
	2.5mm		01
	3.0mm		01
	3.7mm		01
	<b>Recommendably lengths are 15cm, 20cm, 25cm, 30cm, 35cm</b>		
<p><b>3</b></p>  <p><b>FLAT FLAPPER</b> <i>"It is used to open and expose adipose tissue for more precise suctioning"</i></p>	1.5mm		01
	2.0mm		01
	3.0mm		01
	<b>Recommendably lengths are 10cm, 15cm, 20cm</b>		
<p><b>4</b></p>  <p><b>V Tunnel Dissector</b> <i>"Slight upward curve to tip in order to inhibit direct forceful contact, comfortable handle with orientation marking."</i></p>	1.5mm		01
	2.0mm		01
	2.5mm		01
	3.0mm		01
	3.7mm		01
	<b>Recommendably lengths are 7cm, 10cm, 15cm, 20cm, 25cm, 30cm, 35cm</b>		
<p><b>5</b></p>  <p><b>W Dissector Cannulas</b> <i>"Small areas of fibrous tissue, less aggressive and finer than V Dissector."</i></p>	1.5mm		01
	2.0mm		01
	<b>Recommendably lengths are 10cm, 15cm</b>		



**FILTRATOR**  
LIPOSUCTION

	Simple Luer Lock	Super Luer Lock	Tommy	Handle Fixed	Adopter Luer Lock
	46.175120	46.175220			46.181520
	46.176120	46.176220	46.176320	46.176420	46.176520
	46.176125	46.176225	46.176325	46.176425	46.176525
	46.176130	46.176230	46.176330	46.176430	46.176530
	46.176137	46.176237	46.176337	46.176430	46.176535
	46.177115				
		46.177220	46.177320	46.177420	46.177520
		46.177230	46.177330	46.177430	46.177530
	46.178115	46.178215	40.178315	46.178415	
	46.178120	46.178220	46.178320	46.178420	46.178520
	46.178125	46.178225	46.178325	46.178425	46.178525
		46.178230	46.178330	46.178430	46.178530
			46.178335	46.178435	46.178535
	46.179115	46.179215			
	46.179120	46.17220			

