


# Nostril Retainer



- Sterilization Method: Autoclaving, Boiling

## Features

1. This product is suitable for patients with nasal foramens of different sizes or those with a thick columella nasi wall.
2. 6 sizes are available for both left and right parts, so that the proper size can be selected for each patient.
3. Made of soft, less irritative silicone.


## For Right Side

Product No.	Size a×b×c×d	Packing
#8263	C 5.5×9.5×11×13	1pc/box
#8264	D 6×10.5×12.5×15.5	1pc/box
#8265	E 7×11×14×17	1pc/box
#8266	F 8×12×15.5×19.5	1pc/box
#8267	G 9×13.5×17×21.5	1pc/box
#8268	H 10×14.5×18.5×23.5	1pc/box

## For Left Side

Product No.	Size a×b×c×d	Packing
#8283	C 5.5×9.5×11×13	1pc/box
#8284	D 6×10.5×12.5×15.5	1pc/box
#8285	E 7×11×14×17	1pc/box
#8286	F 8×12×15.5×19.5	1pc/box
#8287	G 9×13.5×17×21.5	1pc/box
#8288	H 10×14.5×18.5×23.5	1pc/box

## Mold Guide



Product No.	Product Name	Packing
#8300	Mold Guide	1set/box

\* Mold Guide for proper size selection

## Precautions for use

- Do not damage the product while picking with forceps. Product may be torn from the damaged part.
- If damaged, please immediately switch to a new product.

The product here prohibits the following acts.



Do not reuse



Refrain from applying iodine antiseptics



Refrain from exposing Ultraviolet Ray for long time