


# Nostril Retainer

Nostril Retainer



- Sterilization Method: Autoclaving, Boiling

## Features

1. This product prevents the deformation of the nasal foramen after the surgery for cleft lip or palate, rhinoplasty, or augmentation rhinoplasty without imposing a burden on the patient.
2. It is available in 13 sizes each for the standard and long types so that it can be used in a size appropriate to each patient.
3. It is made of soft, less irritative silicone.


### Standard Type

Product No.	Size W×L×t(mm)	Packing
#8101	16×8×7	1pc/box
#8102	17×9×8	1pc/box
#8103	18×10×8	1pc/box
#8104	19×10×9	1pc/box
#8105	20×10×10	1pc/box
#8106	21×12×11	1pc/box
#8107	22×12×11	1pc/box
#8108	24×13×13	1pc/box
#8109	25×14×14	1pc/box
#8110	26×15×15	1pc/box
#8111	27×16×16	1pc/box
#8112	29×17×17	1pc/box
#8113	30×17×18	1pc/box

### Long Type

Product No.	Size W×L×t(mm)	Packing
#8131	16×14×7	1pc/box
#8132	17×15×8	1pc/box
#8133	18×16×8	1pc/box
#8134	19×16×9	1pc/box
#8135	20×16×10	1pc/box
#8136	21×18×11	1pc/box
#8137	22×18×11	1pc/box
#8138	24×19×13	1pc/box
#8139	25×20×14	1pc/box
#8140	26×21×15	1pc/box
#8141	27×22×16	1pc/box
#8142	29×23×17	1pc/box
#8143	30×23×18	1pc/box

## Mold Guide



Product No.	Product Name	Packing
#8100	Mold Guide for Standard Type	1set/box
#8150	Mold Guide for Long Type	1set/box

\*Mold Guide for proper size selection.

## Precautions for use

- Do not damage the product while picking with forceps. Product may be torn from the damaged part.
- If damaged, please immediately switch to a new product.

The product here prohibits the following acts.



Do not reuse



Refrain from applying iodine antiseptics



Refrain from exposing Ultraviolet Ray for long time