



TECH 3000 INFILTRATION PUMP



Specifications :-

- Speed: 1 to 300 rpm, Reversible/Bi Directional
- Max Flow Rate : 1140ml/min
- Control : Pneumatic Foot Switch
- Speed resolution: 1 rpm
- Speed control: Rotary encoded switch
- Display: 3-digit LED displays current speed
- Power supply: AC 90 - 260V 50/60 Hz
- Dimensions (L×W×H): 285×207×180 (mm)
- Max Flow Rate : 1140ml/min
- Weight 4.1 Kg

Includes :
1 Reusable Silicone Tubing or 5 Disposable Single Spike Infiltration Tubing
1 Pneumatic Foot Switch, 1 Infiltration Handle , 1 Infiltration Cannula



Reusable Silicone Tubing



Standard Infusion Handle



Single Spike Disposable Infusion Tubing
(Sterile)

REGD. OFFICE
TECHNOMED (INDIA) PRIVATE LIMITED

Regd. Off.: H-107, Shivaji Park, Punjabi Bagh, New Delhi-110026

Tel.: +91 11 41420555

Email : technomedindia@outlook.com