

# BiLumix

SHADOWLESS HEADLAMP

Exclusively Distributed by:



[www.technomed.net](http://www.technomed.net)

+91 85869 52121

[technomedindiapvtltd](https://www.instagram.com/technomedindiapvtltd)



**LOUPES**  
1.5x, 2.5x, 3.5x

**BILUMIX HEADLAMPS AND LOUPES**

Be bold. Illuminate your vision.





## Weightless

### Light and Sturdy

No pressure or added weight on the nose. Contoured to fit snugly on the head for a comfortable wear all day long.



## Shadowless

### Dual Lights

Two LED light sources deliver a strong spotlight in angles that eliminate any shadow from a first person POV.



## Wireless

### Two Rechargeable Batteries

An enhanced power management system along with two rechargeable batteries allow you to work continuously for up to 12 hours.



## SPECIFICATIONS

Weight : 140g  
Brightness : Max 100,000 lux

## BATTERY LIFE :

Max Intensity : About 3 hours (Continuously)  
Mid to Low Intensity : 5 to 8 hours (Continuously)







## INCLUDED ACCESSORIES

BiLumix Headlamps comes with  
2 Rechargeable Batteries,  
1 Battery Cradle,  
1 Attached Light,  
Source Filters (Clear) and  
1 Travel Bag



## ACCESSORIES :

1.5 X Loupe  
Headlamp Loupe Specification

Magnification: 1.5x  
Working distance: 5~100cm  
Working field size: 45cm  
Weight: 13g





## ACCESSORIES :

2.5 X Loupes

Headlamp Loupe Specification

Magnification: 2.5x

Working distance: 20 ~ 100cm (Adjustable)

Weight: 45g





## ACCESSORIES :

3.5 X Loupe

Loupe Specification

Magnification: 3.5x

Working Distance: Adjustable (20 ~ 60cm)

Weight: 51g



Exclusively Distributed by:



www.technomed.net

+91 85869 52121

technomedindiapvtltd

GENERATION 2.0

**BiLumix**  
SHADOWLESS HEADLAMP  
Visit [bilumix.com](http://bilumix.com)

**BRIGHT CHOICE**  
For any procedure to illuminate an area via a unique, dual light source that eliminates any shadow

**WEIGHTLESS**  
Light & Sturdy

**SHADOWLESS**  
Dual lights

**WIRELESS**  
2 Rechargeable Batteries

**BiLumix**  
SHADOWLESS HEADLAMP

**Ergonomic Design** Adjustable headlamp length and front pads allows users to customize their Bilumix for a perfect fit.

**Optimized Comfort** Newly designed padded front and back provides a more comfortable grip and snug fit.

**One-touch On-Off**  
Additional on and off switch on back to help you control your light more efficiently during your procedures.

**Power Indicator**  
Easy way to check battery status.

**Easy Control**  
Quickly change brightness of LED with dial control.

## Technomed India Pvt. Ltd.

Regd. Office : H-107, Shivaji Park, Punjabi Bagh, New Delhi-110026

+91 8586952121 | Office : 011 - 41420555

E-mail : [info@technomed.net](mailto:info@technomed.net) | [technomedindia@outlook.com](mailto:technomedindia@outlook.com)

Website : [www.technomed.net](http://www.technomed.net) | [www.pressuregarment.com](http://www.pressuregarment.com)

**BiLumix**  
SHADOWLESS HEADLAMP

41 University Drive Suite 400,

Newtown, PA 18940

United States

Email:[sales@bilumix.com](mailto:sales@bilumix.com)