

BiLumix

SHADOWLESS HEADLAMP

Exclusively Distributed by:



www.technomed.net

www.shoptechnomed.com

+91 85869 52121

[technomedindia Pvt Ltd](#)

BILUMIX HEADLAMPS AND LOUPES



GEN 2

LOUPES
1.5x, 2.5x, 4x



GEN 2S

LOUPES
1.5x, 2.5x, 4x

Be bold. Illuminate your vision.



BiLumix

Illuminate Your Vision

For any procedure to illuminate an area via a unique, dual light source that eliminates any shadow.



NEW
2S

- Adjustable Light Angle
- Focused Light


- Shadowless
- Weightless
- Wireless
- Ergonomic Design



Weightless

Light and Sturdy

No pressure or added weight on the nose. Contoured to fit snugly on the head for a comfortable wear all day long.



Shadowless

Dual Lights

Two LED light sources deliver a strong spotlight in angles that eliminate any shadow from a first person POV.



Wireless

Two Rechargeable Batteries

An enhanced power management system along with two rechargeable batteries allow you to work continuously for up to 12 hours.



SPECIFICATIONS

Weight : 140g
Brightness : Max 100,000 lux

BATTERY LIFE :

Max Intensity : About 3 hours (Continuously)
Mid to Low Intensity : 5 to 8 hours (Continuously)





INCLUDED ACCESSORIES

BiLumix Headlamps comes with
2 Rechargeable Batteries,
1 Battery Cradle,
1 Attached Light,
Source Filters (Clear) and
1 Travel Bag

ACCESSORIES :

1.5 X Loupe
Headlamp Loupe Specification

Magnification: 1.5x
Working distance: 5~100cm
Working field size: 45cm
Weight: 13g





ACCESSORIES :

2.5 X Loupes

Headlamp Loupe Specification

Magnification: 2.5x

Working distance: 20 ~ 100cm (Adjustable)

Weight: 45g





ACCESSORIES :

4 X Loupe

Loupe Specification

Magnification: 4x

Working Distance: Adjustable (20 ~ 60cm)

Weight: 51g

Exclusively Distributed by:



www.technomed.net

www.shoptechnomed.com

+91 85869 52121



technomedindiapvtltd



GENERATION 2.0

BiLumix
SHADOWLESS HEADLAMP

BRIGHT CHOICE
For any procedure to illuminate an area via a unique, dual light source that eliminates any shadow

WEIGHTLESS
Light & Sturdy

SHADOWLESS
Dual Lights

WIRELESS
2 Rechargeable Batteries

BiLumix
SHADOWLESS HEADLAMP

Ergonomic Design Adjustable headlamp length and front pads allows users to customize their Bilumix for a perfect fit.

Optimized Comfort Newly designed padded front and back provides a more comfortable grip and snug fit.



One-touch On-Off
Additional on and off switch on back to help you control your light more efficiently during your procedures.

Power Indicator
Easy way to check battery status.

Easy Control
Quickly change brightness of LED with dial control.

Technomed India Pvt. Ltd.

Regd. Office : C- 19, DSIIDC Complex, Kirti Nagar, New Delhi - 110015

+91 8586952121 | Office : 011 - 41420555

E-mail : info@technomed.net | technomedindia@outlook.com

Website : www.technomed.net | www.pressuregarment.com | www.shoptechnomed.com

BiLumix
SHADOWLESS HEADLAMP

41 University Drive Suite 400,

Newtown, PA 18940

United States

Email:sales@bilumix.com