

Ultrasonic Bone Surgery , Imaging & Endoscopy

# CARING PEOPLE, IMPROVING LIVES.

Medical Catalog

# Clinical kits

## Rhinoplasty tips «Essential»



Reference	Description
F87681	Rhinoplasty kit "Essential" non-sterile 6 tips - 5 times reusable

*Metallic support included*



Reference	F87686	F87687	F87679	F87677	F87678	F87680
Length (mm)	33.7	33.7	34	30	30	27
Offset (mm)	11	11	11	10	10	17
Width (mm)	3.5	3.5	2.6	5	5	•
Thickness (mm)	1.3	0.9	0.5	0.6	0.6	0.5

# Clinical kits

## Rhinoplasty tips «Expert»



Reference	Description
F87689	Rhinoplasty kit "Expert" non-sterile 9 tips - 5 times reusable

*Metallic support included*



Reference	F87688	F87686	F87687	F87677	F87678
Length (mm)	28	33.7	33.7	30	30
Offset (mm)	12.2	11	11	10	10
Width (mm)	3.5	3.5	3.5	5	5
Thickness (mm)	0.55	1.3	0.9	0.6	0.6



Reference	F87683	F87682	F87679	F87680
Length (mm)	31.5	31.5	34	27
Offset (mm)	9.5	9.5	11	17
Width (mm)	2.6	2.6	2.6	•
Thickness (mm)	0.5	0.5	0.5	0.5

# Clinical kits

## Rhinoplasty tips «Preservation»



**RHL4L\*\***      **RHL4R\*\***      **RHL5\*\***      **RHL6\*\***      **BS6 XXL**

Reference	F87692	F87691	F87684	F87693	F87604
Length (mm)	91.7	91.7	59.6	101.4	103
Offset (mm)	9.7	9.7	14.6	9.7	12.2
Width (mm)	2.6	2.6	5.1	5	3.5
Thickness (mm)	0.4	0.4	0.6	0.6	0.55

# Clinical kits

## CMF tips «Essential»



Reference	Description
F57803	CMF kit non-sterile 6 tips - 5 times reusable

*Metallic support included*



**BS1L**      **BS2L XL**      **BS2R XL**      **BS1RD**      **SL1**      **BS4**

Reference	F87612	F87605	F87606	F87608	F87618	F87615
Length (mm)	31	39.5	39.5	38.9	33	31
Offset (mm)	22	30	30	12.3	11	12.6
Width (mm)	4	4	4	5	2.6	4
Thickness (mm)	0.6	0.8	0.8	0.6	0.3	0.7

\*\*Please contact us to find out the availability of the product in your country

# All tips

## Osteotomy tips



**BS1S**      **BS1RD**      **BS1L**      **BS2L**      **BS2R**

Reference	F87611	F87608	F87612	F87613	F87614
Length (mm)	28.2	38.9	31	28	28
Offset (mm)	15	12.3	22	9.35	9.35
Width (mm)	4	5	4	4	4
Thickness (mm)	0.6	0.6	0.6	0.65	0.65



**BS2L XL**      **BS2R XL**      **BS5**      **RHS3L**      **RHS3R**

Reference	F87605	F87606	F87616	F87677	F87678
Length (mm)	39.5	39.5	33.8	30	30
Offset (mm)	30	30	11.1	10	10
Width (mm)	4	4	2.6	5	5
Thickness (mm)	0.8	0.8	0.5	0.6	0.6



**BS1 XXL**      **RHL4L**      **RHL4R**      **RHL5**      **RHL6**

Reference	F87602	F87692	F87691	F87684	F87693
Length (mm)	104	91.7	91.7	59.6	101.4
Offset (mm)	15	9.7	9.7	14.6	9.7
Width (mm)	4	2.6	2.6	5.1	5
Thickness (mm)	0.6	0.4	0.4	0.6	0.6



**RHS4L**      **RHS4R**      **RHS5**      **NINJA**      **LC1 XL**      **SL1**      **CS1**

Reference	F87683	F87682	F87679	F87628	F87607	F87618	F87629
Length (mm)	31.5	31.5	34	26.4	32.1	33	29.8
Offset (mm)	9.5	9.5	11	15.2	21.5	11	17.8
Width (mm)	2.6	2.6	2.6	2.6	1.2	2.6	4
Thickness (mm)	0.5	0.5	0.5	0.6	0.7	0.3	0.5

# All tips

## Osteoplasty tips



	BS4	BS6	RHS1	RHS2Hb	RHS2Fb
<b>Reference</b>	F87615	F87617	F87688	F87686	F87687
Length (mm)	31	28	28	33.7	33.7
Offset (mm)	12.6	11.9	12.2	11	11
Width (mm)	4	3.5	3.5	3.5	3.5
Thickness (mm)	0.7	0.55	0.55	1.3	0.9



BS6 XXL	
<b>Reference</b>	F87604
Length (mm)	103
Offset (mm)	12.2
Width (mm)	3.5
Thickness (mm)	0.55

# All tips

## Bone remodeling / Drilling tips



	CE3	DB1	DB2	RHS6	SL2
<b>Reference</b>	F87653	F87610	F87609	F87680	F87619
Length (mm)	29.8	27	30.7	27	28
Offset (mm)	10.2	17	20.1	17	14
Width (mm)	•	•	•	•	•
Thickness (mm)	0.85	0.6	2.3	0.5	1.5

# All tips

## Membrane detachment tips



	SL3	SL4	SL5
<b>Reference</b>	F87620	F87621	F87622
<b>Length (mm)</b>	31	29.1	30.9
<b>Offset (mm)</b>	12	10.2	11.1
<b>Width (mm)</b>	5	4	4
<b>Thickness (mm)</b>	0.4	0.65	0.65



# Tip measuring guide

	Length (a)	Offset (b)	Width (c)	Thickness (d)
Osteotomy				
Osteoplasty				
Bone remodeling / Drilling				
Membrane detachment				



 **SATELEC**

A company of ACTEON® Group

**17 Av. Gustave Eiffel,  
ZI Phare  
33700 Mérignac  
France**

Follow-us

